

Anti-Git2 Rabbit Polyclonal Antibody For Wb, Elisa - Q14161

Item Number: CM0012840

Introduction

Rabbit polyclonal antibody against human GIT2 (ARF GTPase-activating protein GIT2), validated for Western blot and ELISA. Ideal for studying GPCR signaling and cytoskeleton dynamics. - Q14161

[Learn More](#)

Parameter	Value
Product Name	GIT2 Rabbit pAb
Remarks/Alias	PKL; CAT2; CAT-2; GIT2
Species	Human
Gene ID (Human)	9815
Gene ID	9815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-510 of human GIT2 (NP_476511.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14161
Protein Weight	85kDa
Shipping	Ice bag