


Anti-Gitr Ligand/Tnfsf18 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ung2

Item Number: CM0027827



Introduction

Rabbit polyclonal antibody against human GITR Ligand (TNFSF18), also known as AITRL/GITRL. Validated in WB and ELISA, cross-reacts with mouse and rat. For T-cell activation and immunology studies.

[Learn More](#)

Parameter	Value
Product Name	GitR Ligand/TNFSF18 Rabbit pAb
Remarks/Alias	TL6; AITRL; GITRL; TNLG2A; hGITRL; GITR Ligand/TNFSF18
Species	Human
GeneID	8995
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-177 of human GITR Ligand/TNFSF18 (NP_005083.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNG2
Protein Weight	20kDa
Shipping	Ice bag