

Anti-Gitr/Tnfrsf18 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y5U5

Item Number: CM0016500



Introduction

Rabbit polyclonal antibody raised against human GITR/TNFRSF18 (aa 26-255). Detects human and mouse GITR via WB and ELISA, suitable for studying T-cell activation, NF-kB signaling, and immune checkpoint biology. SWISS: Q9Y5U5.

[Learn More](#)

Parameter	Value
Product Name	GITR/TNFRSF18 Rabbit pAb
Remarks/Alias	AITR; GITR; CD357; GITR-D; ENERGEN; GITR/TNFRSF18
Species	Human
GeneID (Human)	8784
GeneID	8784
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-255 of human GITR/TNFRSF18 (NP_683699.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9Y5U5
Protein Weight	26kDa
Shipping	Ice bag