

Anti-Gja4 Polyclonal Antibody For Wb, Elisa - P35212

Item Number: CM0023301

Introduction

Rabbit polyclonal antibody targeting human GJA4 (Connexin-37), validated for Western blot and ELISA, with cross-reactivity to rat. Ideal for gap junction and intercellular communication research.

[Learn More](#)

Parameter	Value
Product Name	GJA4 Rabbit pAb
Remarks/Alias	CX37; GJA4
Species	Human
GeneID (human)	2701
GeneID	2701
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 229-333 of human GJA4 (NP_002051.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P35212
Protein Weight	37kDa
Shipping	Ice bag