

Anti-Gja5 Polyclonal Antibody For Wb, Elisa - P36382

Item Number: CM0029216

Introduction

Rabbit polyclonal antibody against GJA5 (Connexin-40, CX40) validated for WB and ELISA. Detects human, mouse, and rat GJA5. Ideal for gap junction and cardiac research.

[Learn More](#)

Parameter	Value
Product Name	GJA5 Rabbit pAb
Remarks/Alias	CX40; ATFB11; GJA5
Species	Human
GeneID	2702
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 229-358 of human GJA5 (NP_859054.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36382
Protein Weight	40kDa
Shipping	Ice bag