

Anti-Gja8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P48165

Item Number: CM0023153

Introduction

Rabbit polyclonal antibody targeting human GJA8 (Connexin-50). Suitable for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for eye lens gap junction and cataract research.

[Learn More](#)

Parameter	Value
Product Name	GJA8 Rabbit pAb
Remarks/Alias	CAE; CAE1; CX50; CZP1; MP70; CTRCT1; GJA8
Species	Human
Gene ID (Human)	2703
Gene ID	2703
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48165
Protein Weight	48kDa
Shipping	Ice bag