


Anti-Gjb1 Rabbit Polyclonal Antibody For Wb And Elisa - P08034

Item Number: CM0017407



Introduction

Rabbit polyclonal antibody against human GJB1 (Connexin-32) for Western blotting and ELISA. Recognizes human, mouse, and rat GJB1. Ideal for gap junction, Charcot-Marie-Tooth disease, and neuropathy research.

[Learn More](#)

Parameter	Value
Product Name	GJB1 Rabbit pAb
Remarks/Alias	CMTX; CX32; CMTX1; GJB1
Species	Human
GeneID (Human)	2705
GeneID	2705
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 214-283 of human GJB1 (NP_000157.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08034
Protein Weight	32kDa
Shipping	Ice bag