

Anti-Gjb3 Polyclonal Antibody For Wb, Elisa - 075712

Item Number: CM0026459

Introduction

Rabbit polyclonal antibody targeting GJB3 (Connexin-31). Validated for Western blot and ELISA, with reactivity across human, mouse, and rat. Suitable for hearing loss and skin disease research.

[Learn More](#)

Parameter	Value
Product Name	GJB3 Rabbit pAb
Remarks/Alias	EKV; CX31; DFNA2; EKV1; DFNA2B; GJB3
Species	Human
GeneID (Human)	2707
GeneID	2707
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human GJB3 (NP_076872.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075712
Protein Weight	31kDa
Shipping	Ice bag