


Anti-Gjd2 Rabbit Pab - Q9Ukl4

Item Number: CM0030873



Introduction

GJD2 Rabbit pAb detects gap junction delta-2 (connexin-36) in human, mouse, rat via WB, ELISA. Target expressed in neuronal gap junctions, ~36kDa. Suitable for neuroscience and cell communication studies.

[Learn More](#)

Parameter	Value
Product Name	GJD2 Rabbit pAb
Remarks/Alias	CX36; GJA9; GJD2
Species	Human
Gene ID (Human)	57369
Gene ID	57369
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-197 of human GJD2 (NP_065711.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKL4
Protein Weight	36kDa
Shipping	Ice bag