

Anti-Gk2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q14410

Item Number: CM0022512

Introduction

Rabbit polyclonal antibody against human GK2 (Glycerol kinase 2, testis-specific). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Useful for male fertility and spermatogenesis research, including mitochondrial sheath formation.

[Learn More](#)

Parameter	Value
Product Name	GK2 Rabbit pAb
Remarks/Alias	GKP2; GKTA; GK2
Species	Human
GeneID (Human)	2712
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14410
Protein Weight	61kDa
Shipping	Ice bag