

Anti-GLB1L3 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nci6

Item Number: CM0015957

Introduction

This rabbit polyclonal antibody against human GLB1L3 (beta-galactosidase-1-like protein 3) is validated for Western blot and ELISA, and cross-reacts with mouse and rat. An essential tool for investigating glycosidase activity and ganglioside metabolism.

[Learn More](#)

Parameter	Value
Product Name	GLB1L3 Rabbit pAb
Remarks/Alias	GLB1L3
Species	Human
GeneID (Human)	112937
GeneID	112937
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-290 of human GLB1L3 (NP_001073876.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NCI6
Protein Weight	75kDa
Shipping	Ice bag