

Anti-Gldc Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P23378

Item Number: CM0029178



Introduction

Rabbit polyclonal antibody against human GLDC, validated for WB, IF/ICC, IP, and ELISA assays. Cross-reacts with human, mouse, and rat samples. Ideal for glycine metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GLDC Rabbit pAb
Remarks/Alias	GCE; GCSP; HYGN1; GLDC
Species	Human
GeneID (Human)	2731
GeneID	2731
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-290 of human GLDC (NP_000161.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23378
Protein Weight	113kDa
Shipping	Ice bag