

Anti-Gldc Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P23378

Item Number: CM0006902



Introduction

Recombinant rabbit monoclonal antibody targeting human GLDC (GCE/GCSP), the mitochondrial P protein of the glycine cleavage system. Suitable for WB, IHC-P, ELISA and cross-reacts with mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	GLDC Rabbit mAb
Remarks/Alias	GCE; GCSP; HYGNI
Species	Human
Gene ID (Human)	2731
Gene ID	2731
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 580-749 of human GLDC (NP_000161.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23378
Protein Weight	113kDa
Shipping	Ice bag