

Anti-Gldc Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - P23378

Item Number: CM0019085



Introduction

Rabbit polyclonal antibody against human GLDC, suitable for WB, IF, IHC-P, and ELISA. Validated in human, mouse, rat. Targets mitochondrial glycine decarboxylase.

[Learn More](#)

Parameter	Value
Product Name	GLDC Rabbit pAb
Remarks/Alias	GLDC; GCE; GCSP; HYGN1; glycine decarboxylase
Species	Human
GeneID (Human)	2731
GeneID	2731
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 580-749 of human GLDC (NP_000161.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23378
Protein Weight	112kDa
Shipping	Ice bag