

Anti-Gle1 Polyclonal Antibody For Wb, Elisa - Q53Gs7

Item Number: CM0015325



Introduction

Rabbit anti-human GLE1 polyclonal antibody validated for Western blot and ELISA. Enables investigation of mRNA export mechanisms and nuclear pore complex biology. Essential for nucleoporin research.

[Learn More](#)

Parameter	Value
Product Name	GLE1 Rabbit pAb
Remarks/Alias	LCCS; CAAHC; CAAHD; GLE1L; LCCS1; hGLE1; GLE1
Species	Human
GeneID (Human)	2733
GeneID	2733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GLE1 (NP_001003722.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q53GS7
Protein Weight	80kDa
Shipping	Ice bag