

Anti-Gli1 Rabbit Polyclonal Antibody For Wb And Elisa - P08151

Item Number: CM0031055



Introduction

High-quality anti-GLI1 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Detects human, mouse, and rat GLI1, a key Hedgehog signaling transcription factor implicated in development and cancer.

[Learn More](#)

Parameter	Value
Product Name	GLI1 Rabbit pAb
Remarks/Alias	GLI; PPD1; PAPA8; GLI1
Species	Human
Gene ID	2735
Immunogen	Recombinant protein of human GLI1
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08151
Protein Weight	118kDa
Shipping	Ice bag