

Anti-Gli1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08151

Item Number: CM0016348



Introduction

Rabbit polyclonal antibody against human GLI1, validated for WB, IHC-P, IF/ICC, and ELISA. GLI1 is a key Hedgehog pathway transcription activator implicated in development and cancer. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	GLI1 Rabbit pAb
Remarks/Alias	GLI; PPD1; PAPA8; GLI1
Species	Human
Gene ID (Human)	2735
Gene ID	2735
Immunogen	Synthetic peptide corresponding to amino acids 350-450 of human GLI1 (NP_001153517.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08151
Protein Weight	118kDa
Shipping	Ice bag