

Anti-Gli2 Monoclonal Antibody For If/Icc, If-P, Ihc-P, Elisa - P10070

Item Number: CM0006773

Introduction

Recombinant rabbit monoclonal antibody against GLI2, validated for IF/ICC, IF-P, IHC-P, and ELISA in human, mouse, and rat. Suitable for hedgehog signaling pathway research and HTLV-1 studies.

[Learn More](#)

Parameter	Value
Product Name	GLI2 Rabbit mAb
Remarks/Alias	CJS; HPE9; PHS2; THP1; THP2; GLI2
Species	Human
GeneID (human)	2736
GeneID	2736
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1100-1200 of human GLI2. (NP_005261.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10070
Shipping	Ice bag