

Anti-Gli2 Polyclonal Antibody For Wb, Elisa - P10070

Item Number: CM0026024

Introduction

Rabbit polyclonal antibody against GLI2, suitable for WB and ELISA, with cross-reactivity to mouse and rat. Targets transcription factor GLI2 involved in hedgehog signaling and development. For research use only.

[Learn More](#)

Parameter	Value
Product Name	GLI2 Rabbit pAb
Remarks/Alias	GLI2; CJS; HPE9; PHS2; THP1; THP2; Gli2
Species	Mouse
GeneID (Human)	2736
GeneID	14633
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of mus musculus GLI2 (NP_001074594.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P10070
Protein Weight	168kDa
Shipping	Ice bag