

Anti-Gli2 Polyclonal Antibody For Wb, Elisa - P10070

Item Number: CM0024629

Introduction

Rabbit polyclonal antibody targeting GLI2 (Gene ID 2736, SWISS P10070) validated for WB and ELISA. Reacts with human, mouse, rat. Suitable for smoothened signaling and transcription factor research.

[Learn More](#)

Parameter	Value
Product Name	GLI2 Rabbit pAb
Remarks/Alias	CJS; HPE9; PHS2; THP1; THP2; GLI2
Species	Human
Gene ID	2736
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1100-1200 of human GLI2 (NP_005261.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10070
Protein Weight	168kDa
Shipping	Ice bag