


Anti-Gli3 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P10071

Item Number: CM0016296



Introduction

Rabbit polyclonal antibody against human GLI3, validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for studying hedgehog signaling and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	GLI3 Rabbit pAb
Remarks/Alias	PHS; ACLS; GCPS; PAPA; PAPB; PAP-A; PAPA1; PPDIV; GLI3FL; GLI3-190; GLI3
Species	Human
GeneID (Human)	2737
GeneID	2737
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 681-788 of human GLI3 (NP_000159.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10071
Protein Weight	170kDa
Shipping	Ice bag