

Anti-Glis3 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nea6

Item Number: CM0017855

Introduction

Anti-GLIS3 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat GLIS3. Ideal for transcription and diabetes research. - Q8NEA6.

[Learn More](#)

Parameter	Specification
Product Name	GLIS3 Rabbit pAb
Remarks/Alias	NDH; ZNF515; GLIS3
Species	Human
GeneID (Human)	169792
GeneID	169792
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 606-775 of human GLIS3 (NP_689842.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NEA6
Protein Weight	84kDa
Shipping	Ice bag