

Anti-Glo1 Monoclonal Antibody For Wb, Elisa - Q04760

Item Number: CM0010507



Introduction

Rabbit monoclonal antibody against human GLO1 (Lactoylglutathione lyase), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying glyoxalase detoxification and NF-κB signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	GLO1 Rabbit mAb
Remarks/Alias	GLYI; GLOD1; HEL-S-74; GLO1
Species	Human
Gene ID (Human)	2739
Gene ID	2739
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GLO1 (Q04760).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04760
Protein Weight	21kDa
Shipping	Ice bag