

Anti-Glod4 Polyclonal Antibody For Wb, Elisa - Q9Hc38

Item Number: CM0027538



Introduction

Rabbit polyclonal antibody targeting human GLOD4 (glyoxalase domain-containing protein 4). Validated for Western blot and ELISA, with cross-reactivity to mouse. Suitable for mitochondrial research and glyoxalase studies.

[Learn More](#)

Parameter	Value
Product Name	GLOD4 Rabbit pAb
Remarks/Alias	HC6; HC71; CGI-150; C17orf25; GLOD4
Species	Human
GeneID (Human)	51031
GeneID	51031
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 179-298 of human GLOD4 (NP_057164.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HC38
Protein Weight	35kDa
Shipping	Ice bag