

Anti-Glomulin (Glmn) Polyclonal Antibody For Wb, Elisa - Q92990

Item Number: CM0017295



Introduction

Rabbit polyclonal antibody against Glomulin (GLMN), validated for WB and ELISA. Cross-reacts with human and mouse. Essential for SCF ubiquitin ligase regulation and vasculature development. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	GLMN Rabbit pAb
Remarks/Alias	FAP; GVM; GLML; FAP48; FAP68; FKBPAP; VMGLOM; GLMN
Species	Human
GeneID (Human)	11146
GeneID	11146
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 405-594 of human GLMN (NP_444504.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92990
Protein Weight	68kDa
Shipping	Ice bag