


Anti-Glp-1 Monoclonal Antibody For If-P, Ihc-P, Elisa - P01275

Item Number: CM0006083



Introduction

High-quality rabbit monoclonal antibody targeting human GLP-1. Validated for immunofluorescence, immunohistochemistry, and ELISA. Cross-reactive with mouse and rat. Supports diabetes, obesity, and metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	GLP-1 Rabbit mAb
Remarks/Alias	GLP1; GLP2; GRPP; GLP-1
Species	Human
GeneID	2641 (Human)
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 70-150 of human GLP-1 (P01275).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P01275
Protein Weight	21kDa
Shipping	Ice bag