


Anti-Glp1R Polyclonal Antibody For Wb And Elisa - P43220

Item Number: CM0017491



Introduction

High-quality rabbit polyclonal antibody targeting human GLP1R (GLP-1 receptor). Validated for WB and ELISA, cross-reacts with human and rat. Suitable for diabetes, obesity, and metabolic disorder research.

[Learn More](#)

Parameter	Specification
Product Name	GLP1R Rabbit pAb
Remarks/Alias	GLP-1; GLP-1R; GLP-1-R; GLP1R
Species	Human
Gene ID	2740 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-116 of human GLP1R (NP_002053.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P43220
Protein Weight	53kDa
Shipping	Ice bag