


Anti-Glp1R Polyclonal Antibody For Wb, If/Icc, Elisa - P43220

Item Number: CM0030131



Introduction

Rabbit polyclonal antibody against human GLP1R, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for glucagon-like peptide-1 receptor signaling and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	GLP1R Rabbit pAb
Remarks/Alias	GLP-1; GLP-1R; GLP-1-R; GLP1R
Species	Human
GeneID (Human)	2740
GeneID	2740
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human GLP1R (NP_002053.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P43220
Protein Weight	53kDa
Shipping	Ice bag