


Anti-Glp2R Monoclonal Antibody For Wb And Elisa - O95838

Item Number: CM0003188



Introduction

Rabbit monoclonal antibody targeting human GLP2R (Glucagon-like peptide 2 receptor) for WB and ELISA. Cross-reacts with mouse and rat. Ideal for GPCR signaling and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	GLP2R Rabbit mAb
Remarks/Alias	GLP2R
Species	Human
Gene ID (Human)	9340
Gene ID	9340
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-553 of human GLP2R (NP_004237.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95838
Protein Weight	63kDa
Shipping	Ice bag