


Anti-Glp2R Polyclonal Antibody For Wb, Ihc-P, Elisa - 095838

Item Number: CM0027724



Introduction

CamelBio provides anti-GLP2R polyclonal antibody, validated for WB, IHC-P, and ELISA. Recognizes human and mouse GLP2R (63 kDa), a GPCR involved in intestinal hormone signaling. Ideal for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	GLP2R Rabbit pAb
Remarks/Alias	GLP2R
Species	Human
Gene ID (Human)	9340
Gene ID	9340
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 460-553 of human GLP2R (NP_004237.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95838
Protein Weight	63kDa
Shipping	Ice bag