

Anti-Glra2 Rabbit Polyclonal Antibody For Wb And Elisa - P23416

Item Number: CM0014074



Introduction

Rabbit polyclonal antibody targeting human GLRA2 (glycine receptor subunit alpha-2). Validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat. Ideal for neuroscience and inhibitory synapse research.

[Learn More](#)

Parameter	Value
Product Name	GLRA2 Rabbit pAb
Remarks/Alias	GLR; MRXSP; GLRA2
Species	Human
GeneID (human)	2742
GeneID	2742
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-420 of human GLRA2 (NP_002054.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23416
Protein Weight	52kDa
Shipping	Ice bag