


Anti-Grx Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35754

Item Number: CM0027202



Introduction

GLRX Rabbit pAb targets human glutaredoxin-1, with cross-reactivity to mouse and rat. Validated for WB, IHC-P, IF/ICC, and ELISA, this polyclonal antibody supports redox biology and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	GLRX Rabbit pAb
Remarks/Alias	GRX; GRX1; GLRX
Species	Human
Gene ID (Human)	2745
Gene ID	2745
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of human GLRX (NP_002055.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35754
Protein Weight	12kDa
Shipping	Ice bag