


Anti-Glrx3 Polyclonal Antibody For Wb, Ihc-P, Elisa - O76003

Item Number: CM0012889



Introduction

Rabbit polyclonal antibody against human GLRX3 (PICOT, GRX4), validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Targets the cytosolic iron-sulfur cluster assembly factor involved in cardiac hypertrophy regulation and hemoglobin maturation.

[Learn More](#)

Parameter	Value
Product Name	GLRX3 Rabbit pAb
Remarks/Alias	GRX3; GRX4; GLRX4; PICOT; TXNL2; TXNL3; GLRX3
Species	Human
GeneID (human)	10539
GeneID	10539
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-335 of human GLRX3 (NP_001186797.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O76003
Protein Weight	37kDa
Shipping	Ice bag