


Anti-Glrx3 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - O76003

Item Number: CM0024836



Introduction

Anti-GLRX3 Rabbit Polyclonal Antibody validated for WB, IP, ELISA in human samples.

Glutaredoxin-3 functions in cytosolic iron-sulfur cluster assembly, cardiac hypertrophy regulation, and hemoglobin maturation. Ideal for cardiovascular and hematology research.

[Learn More](#)

Parameter	Value
Product Name	GLRX3 Rabbit pAb
Remarks/Alias	GRX3, GRX4, GLRX4, PICOT, TXNL2, TXNL3, GLRX3
Species	Human
Gene ID (Human)	10539
Gene ID	10539
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-335 of human GLRX3 (NP_001186797.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	O76003
Protein Weight	37kDa
Shipping	Ice bag