

Anti-Gls2 Polyclonal Antibody For Wb, Elisa - Q9Ui32

Item Number: CM0024510



Introduction

High-quality Anti-GLS2 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, rat. Mitochondrial glutaminase involved in metabolism and tumor suppression. Ideal for IVD development.

[Learn More](#)

Parameter	Value
Product Name	GLS2 Rabbit pAb
Remarks/Alias	GA; GLS; LGA; hLGA; GLS2
Species	Human
GeneID (Human)	27165
GeneID	27165
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 545-602 of human GLS2 (NP_037399.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Ui32
Protein Weight	66kDa
Shipping	Ice bag