


Anti-Glucagon Monoclonal Antibody For If-P, Ihc-P, Elisa - P01275

Item Number: CM0007091



Introduction

Rabbit monoclonal antibody against human Glucagon, validated for IF-P, IHC-P, and ELISA. Reacts with human, mouse, rat. A key tool for studying glucagon signaling in glucose metabolism and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	Glucagon Rabbit mAb
Remarks/Alias	GLP1; GLP2; GRPP; GLP-1; Glucagon
Species	Human
GeneID (Human)	2641
GeneID	2641
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Glucagon (NP_002045.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01275
Protein Weight	21kDa
Shipping	Ice bag