


Anti-Glucagon Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P01275

Item Number: CM0003646



Introduction

High-specificity rabbit monoclonal antibody targeting human Pro-glucagon (Glucagon, GLP-1, GLP-2). Validated for IF-P, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for diabetes and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Glucagon Rabbit mAb
Remarks/Alias	GLP1; GLP2; GRPP; GLP-1; Glucagon
Species	Human
GeneID (Human)	2641
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Glucagon (NP_002045.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P01275
Protein Weight	21kDa
Shipping	Ice bag