


Anti-Glucagon Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01275

Item Number: CM0005402



Introduction

High-quality rabbit monoclonal antibody against human Glucagon (GCG) for WB, IF-P, IHC-P, ELISA. Validated in human, mouse, rat. Ideal for glucose metabolism, diabetes, and gut hormone research.

[Learn More](#)

Parameter	Value
Product Name	Glucagon Rabbit mAb
Remarks/Alias	GLP1, GLP2, GRPP, GLP-1, Glucagon
Species	Human
Gene ID (Human)	2641
Gene ID	2641
Immunogen	Synthetic peptide corresponding to amino acids 21-120 of human Glucagon (NP_002045.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01275
Protein Weight	21kDa
Shipping	Ice bag