


Anti-Glucagon Receptor (Gcgr) Polyclonal Antibody For Wb, If-P, Elisa - P47871

Item Number: CM0017362



Introduction

Rabbit polyclonal antibody against human Glucagon Receptor (GCGR) for WB, IF-P, ELISA. Validated in human, mouse, rat. Suitable for diabetes and glucose metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GCGR Rabbit pAb
Remarks/Alias	GGR; GL-R; MVAH; GCGR
Species	Human
GeneID (Human)	2642
GeneID	2642
Immunogen	Synthetic peptide corresponding to amino acids 50-150 of human GCGR (NP_000151.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47871
Protein Weight	54kDa
Shipping	Ice bag