

Anti-Glucose-6-Phosphate Isomerase (Gpi) Monoclonal Antibody [Kd Validated] For Wb, If/Icc, Elisa - P06744

Item Number: CM0010201



Introduction

KD-validated rabbit monoclonal antibody targeting glucose-6-phosphate isomerase (GPI), also known as AMF/NLK, for WB, IF/ICC, and ELISA. Detects endogenous GPI in human, mouse, and rat samples; ideal for glycolysis and cytokine research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] GPI Rabbit mAb
Remarks/Alias	AMF; NLK; PGI; PHI; GNPI; SA36; SA-36
Species	Human
GeneID (Human)	2821
GeneID	2821
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human glucose-6-phosphate isomerase (NP_000166.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06744
Protein Weight	63 kDa
Shipping	Ice bag