

Anti-Glucose-6-Phosphate Isomerase Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P06744

Item Number: CM0028937



Introduction

Rabbit polyclonal antibody targeting human glucose-6-phosphate isomerase (GPI/AMF/NLK). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for glycolysis, cancer, and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	glucose-6-phosphate isomerase Rabbit pAb
Remarks/Alias	AMF; NLK; PGI; PHI; GNPI; SA36; SA-36; glucose-6-phosphate isomerase
Species	Human
Gene ID (Human)	2821
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human glucose-6-phosphate isomerase (NP_000166.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06744
Protein Weight	63kDa
Shipping	Ice bag