

Anti-Glucosylceramidase Beta (Gba) Polyclonal Antibody For Wb, Elisa - P04062

Item Number: CM0024572



Introduction

Rabbit polyclonal antibody against human glucosylceramidase beta (GBA), suitable for WB and ELISA. Cross-reacts with human and rat. Targets lysosomal lipid degradation enzyme implicated in Gaucher disease.

[Learn More](#)

Parameter	Value
Product Name	Glucosylceramidase beta (GBA) Rabbit pAb
Remarks/Alias	GBA; GCB; GLUC; Glucosylceramidase beta (GBA)
Species	Human
GeneID (human)	2629
GeneID	2629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-250 of human Glucosylceramidase beta (GBA) (NP_000148.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P04062
Protein Weight	60kDa
Shipping	Ice bag