

Anti-Glucosylceramidase Beta (Gba) Rabbit Polyclonal Antibody For Wb And Elisa - P04062

Item Number: CM0026444



Introduction

Rabbit anti-human GBA polyclonal antibody validated for Western Blot and ELISA. Recognizes lysosomal acid glucosylceramidase (60 kDa), a key enzyme in sphingolipid metabolism, linked to Gaucher and Parkinson's disease.

[Learn More](#)

Parameter	Value
Product Name	Glucosylceramidase beta (GBA) Rabbit pAb
Remarks/Alias	GBA; GCB; GLUC; Glucosylceramidase beta (GBA)
Species	Human
GeneID (Human)	2629
GeneID	2629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-250 of human Glucosylceramidase beta (GBA) (NP_000148.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04062
Protein Weight	60kDa
Shipping	Ice bag