


Anti-GLUD1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00367

Item Number: CM0011966



Introduction

Recombinant rabbit monoclonal antibody targeting human GLUD1 (glutamate dehydrogenase 1). Validated in WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for metabolic and cancer biology studies.

[Learn More](#)

Parameter	Value
Product Name	GLUD1 Rabbit mAb
Remarks/Alias	GDH; GDH1; GLUD; hGDH1; GLUD1
Species	Human
GeneID (human)	2746
GeneID	2746
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 459-558 of human GLUD1 (P00367)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00367
Protein Weight	61kDa
Shipping	Ice bag