


Anti-Glud1/2 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P00367/P49448

Item Number: CM0006457



Introduction

Rabbit monoclonal antibody against GLUD1/2 (GDH1/2) for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Validated for research on glutamate metabolism and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	GLUD1/2 Rabbit mAb
Remarks/Alias	GLUD2; GDH2; GLUDP1; GLUD1/2
Species	Human
GeneID (Human)	2746/2747
GeneID	2746/2747
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 319-558 of human GLUD1/2 (NP_036216.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00367/P49448
Protein Weight	61kDa
Shipping	Ice bag