


Anti-Glud2 Rabbit Polyclonal Antibody - P49448

Item Number: CM0029963



Introduction

Anti-GLUD2 Rabbit Polyclonal Antibody for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat GLUD2. Targets mitochondrial glutamate dehydrogenase 2 involved in glutamate recycling during neurotransmission.

[Learn More](#)

Parameter	Value
Product Name	GLUD2 Rabbit pAb
Remarks/Alias	GDH2; GLUDP1; GLUD2
Species	Human
GeneID (Human)	2747
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 319-558 of human GLUD2 (NP_036216.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49448
Protein Weight	61kDa
Shipping	Ice bag