


Anti-GluR1/Gria1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P42261

Item Number: CM0022268



Introduction

Rabbit polyclonal antibody against human GluR1/GRIA1 (Glutamate receptor 1). Validated for WB, IHC-P, ELISA in human, mouse, rat. Recognizes ~102 kDa protein. Ideal for studying excitatory neurotransmission and synaptic plasticity.

[Learn More](#)

Parameter	Value
Product Name	GluR1/GRIA1 Rabbit pAb
Remarks/Alias	GLUH1; GLUR1; GLURA; GluA1; HBGR1; MRD67; MRT76; GluR1/GRIA1
Species	Human
GeneID (Human)	2890
GeneID	2890
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 330-470 of human GluR1/GRIA1 (NP_000818.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42261
Protein Weight	102kDa
Shipping	Ice bag