


Anti-GluR1/Gria1 Rabbit Monoclonal Antibody For Wb, Elisa - P42261

Item Number: CM0002723



Introduction

High-quality rabbit monoclonal antibody targeting human GluR1/GRIA1 (AMPA receptor subunit). Suitable for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for synaptic function and neurotransmission studies.

[Learn More](#)

Parameter	Value
Product Name	GluR1/GRIA1 Rabbit mAb
Remarks/Alias	GLUH1; GLUR1; GLURA; GluA1; HBGR1; MRD67; MRT76; GluR1/GRIA1
Species	Human
Gene ID	2890
Immunogen	Synthetic peptide corresponding to amino acids 800-906 of human GluR1/GRIA1 (P42261)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P42261
Protein Weight	102kDa
Shipping	Ice bag