


Anti-GluR2+GluR3 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P42262 / P42263

Item Number: CM0004075



Introduction

Rabbit monoclonal antibody targeting human GluR2 and GluR3 AMPA receptor subunits. Validated for WB, IP, ELISA with cross-reactivity to mouse and rat. Suitable for neuroscience studies of synaptic plasticity and excitatory transmission.

[Learn More](#)

Parameter	Value
Product Name	GluR2+GluR3 Rabbit mAb
Remarks/Alias	GLUR2; GLURB; GluA2; HBGR2; NEDLIB; gluR-2; gluR-B; GluR-K2; GLUR3; GLURC; GluA3; MRX94; MRXSW; GLUR-C; iGluR3; GLUR-K3; GluR2+GluR3
Species	Human
GeneID (Human)	2891/2892
GeneID	2891/2892
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 784-883 of human GluR2 (P42262)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P42262/P42263
Protein Weight	98kDa
Shipping	Ice bag