


Anti-GluR2/Gria2 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P42262

Item Number: CM0003435



Introduction

Recombinant rabbit monoclonal antibody against human GluR2/GRIA2 (AMPA receptor subunit). Validated for WB, IF-P, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for neuroscience research and synaptic transmission studies.

[Learn More](#)

Parameter	Value
Product Name	GluR2/GRIA2 Rabbit mAb
Remarks/Alias	GLUR2; GLURB; GluA2; HBGR2; NEDLIB; glur-2; glur-B; GluR-K2; GluR2/GRIA2
Species	Human
Gene ID (Human)	2891
Gene ID	2891
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human GluR2/GRIA2 (P42262)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P42262
Protein Weight	99kDa
Shipping	Ice bag